# Levenhuk 5S/7S NG Microscopes

User Manual  
Návod k použití  
Bedienungsanleitung  
Guía del usuario  
Használati útmutató  
Instrukcja obsługi  
Инструкция по эксплуатации


*Radost zaostřit  
Zoom ran und hab Fun!  
Amplie y disfrute  
Kellemes nagyítást!  
Radość przybliżania  
Приближает с удовольствием*

**levenhuk**<sup>®</sup>  
Zoom&Joy

Levenhuk 5S NG

Levenhuk 7S NG

**EN**

- 1 Monocular head (eyepiece tube)
- 2 Revolving nosepiece
- 3 Objectives
- 4 Iris diaphragm (7S NG)
- 5 Mirror
- 6 Coarse adjustment knob
- 7 Fine adjustment knob
- 8 Stage with specimen holders (5S NG)
- 9 Stage with mechanical scale (7S NG)
- 10 Base
- 11 Abbe condenser (7S NG)
- 12 Diaphragm disk (5S NG)

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These microscopes are designed for general biology and bacteriology research in schools, colleges, hospitals, labs, research institutions and universities as well as for clinical testing and teaching demonstrations. By choosing a different objective and eyepiece, different total magnifications can be obtained. The maximum magnifications are 500x (for 5S NG model) and 800x (for 7S NG model).

## Using the microscope

1. Place the specimen slide on the stage and clamp it with clips.
2. Install objectives into the revolving nosepiece, and an eyepiece into the eyepiece tube, respectively.
3. Use a low power objective to observe the object first, then move the object to the center of the bright field. When the object is centered in the field of view, switch to a high power objective to observe the object in detail.
4. When focusing, first adjust the focusing knob forward until a sharp image is obtained. When using the high power objective, turn the focusing knob slowly and carefully, to prevent the objective from contacting the specimen slide to avoid any possible damage.
5. Turn the reflecting mirror to reflect the light to the stage until the bright field of view is visible, then adjust the aperture of iris diaphragm (disc) for a clear image.

## Care and maintenance

- **Never, under any circumstances, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS.**
- Take necessary precautions when using the device with children or others who have not read or who do not fully understand these instructions.
- After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection.
- Do not try to disassemble the device on your own for any reason. For repairs and cleaning of any kind, please contact your local specialized service center.
- Protect the device from sudden impact and excessive mechanical force. Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws.
- Do not touch the optical surfaces with your fingers. To clean the device exterior, use only special cleaning wipes and special optics cleaning tools from Levenhuk. Do not use any corrosive or acetone-based fluids to clean the optics.
- Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush.
- Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity.
- Be careful during your observations, always replace the dust cover after you are finished with observations to protect the device from dust and stains.
- If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope.
- Store the device in a dry, cool place away from hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures.
- When using the microscope, try not to use it near flammable materials or substances (benzene, paper, cardboard, plastic, etc.), as the base may heat up during use, and might become a fire hazard.
- **Seek medical advice immediately if a small part or a battery is swallowed.**

## Specifications

### Levenhuk 5S NG

### Levenhuk 7S NG

	monocular vertical tube	
Optics material	glass	
Magnification, x	40–500	40–800
Eyepiece tube diameter, mm	23.2	
Eyepieces	10x, 12.5x	10x, 20x
Objectives	4x, 10x, 40xs	
Nosepiece	triple	
Stage, mm	110x120, with stage clips	120x120, with measurements
Condenser	diaphragm disk (5 holes)	Abbe condenser NA0.65 with iris diaphragm and 3 filters: Blue, Green, Yellow
Focus	coarse: 50mm; fine: 1.8mm-2.2mm	
Body material	metal	
Illumination	mirror	

The manufacturer reserves the right to make changes to the product range and specifications without prior notice.

## Levenhuk International Lifetime Warranty

All Levenhuk telescopes, microscopes, binoculars and other optical products, except for accessories, carry a **lifetime warranty** against defects in materials and workmanship. Lifetime warranty is a guarantee on the lifetime of the product on the market. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from date of retail purchase. The warranty entitles you to free repair or replacement of the Levenhuk product in any country where a Levenhuk office is located if all warranty conditions are met. For further details please visit our web site: <https://www.levenhuk.com/warranty>  
If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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Levenhuk® is registered trademark of Levenhuk, Inc.
© 2018 Levenhuk, Inc. All rights reserved.
Levenhuk, Inc. 924-D East 124th Ave. Tampa, FL 33612 USA